



surgicalpouch.com



A practical, safe sterilisation pouch for use by sterile service in conjunction with current sized instrument baskets. CE Marked to standards EN11607 1&2, DIN58953-6 Packaging for terminally sterilised medical devices.

Developed with several NHS trusts it acts as a bio barrier protecting sterile surgical instruments through advanced textile technology. Made from 3-layer laminated fabric which protects the contents and will not rip or tear.



Reusable for 200 processes it offers substantial cost savings, is environmentally sustainable and eliminates current surgical wrap which is non-biodegradable. Pouches can remain sterile for up to 12 months using steam sterilisation and 24 months using EO process.

Our products are manufactured in the UK and come complete with 2x Veloclip and security tags. Design features include an open end to allow for thorough cleaning, a GSI Barcode for tracking and monitoring usage and a fixed ID pocket for contents information.

The pouches are compliant to the Health Technical Memorandum 01-01, meaning they keep instruments moist following operations to ensure compliant cleaning processes and can withstand the autoclave to effectively sterilise the equipment. We can produce bespoke orders of any size, colour.

CONTACT US FOR MORE INFORMATION:

 07872 383 183

 info@surgicalpouch.com

 www.surgicalpouch.com

